

## CV of Shanrong Zhang

Senior Manager  
Pre-Clinical Imaging Core  
Center for Biometric Analysis  
The Jackson Laboratory

July 14<sup>th</sup>, 2021

### 1. PERSONAL DATA

<b>Name</b>	Shanrong Zhang, Ph.D.
<b>Home Address</b>	714 State Highway 3 Bar Harbor, ME 04609 (214) 601-2371
<b>Work Address</b>	Pre-Clinical Imaging Center for Biometric Analysis The Jackson Laboratory 600 Main Street Bar Harbor, ME 04609 <a href="mailto:Shanrong.Zhang@jax.org">Shanrong.Zhang@jax.org</a> (207) 288-6974

### 2. EDUCATION

<b>Sep. 1<sup>st</sup> 1990 – Oct. 15<sup>th</sup>, 1995</b>	PhD – Changchun Institute of Applied Chemistry, Chinese Academy of Sciences, Changchun, P.R. China Advisors: Profs. Jiazuan Ni and Fengkui Pei “NMR Studies of Lanthanide Coordination Chemistry”
<b>Sep. 1<sup>st</sup>, 1983 – Jul. 31<sup>st</sup>, 1987</b>	Bachelor of Science, Chemistry Hanzhong Normal College

### 3. EXPERTISE

- Small animal MRI techniques, including the experimental design, the hands-on data acquisition and the image data analysis.
- Expertise in Varian/Agilent VnmrJ 4.0 imaging pulse sequence programming.
- Expertise in Bruker ParaVision 6.0.1 imaging pulse sequence programming.
- Excellent hands-on skills (15+ years) in handling rodents (mouse and rat), including *i.v.*, *s.c.*, *i.m.*, and *i.p.* injection, anesthesia, respiratory / ECG gating, and etc.
- Strong background in MRI contrast agents, from design, development, to in vivo application.
- Excellent skills in nuclear magnetic resonance spectroscopy (NMR), including multi-nuclear (<sup>1</sup>H, <sup>17</sup>O, <sup>31</sup>P, <sup>19</sup>F, <sup>23</sup>Na, etc) and multi-dimensional (2D) NMR.
- Excellent skills in using neuroimaging processing software, such as FSL, ANTs, Brain Connectivity Network, ITK-SNAP, ImageJ, DSISTudio.

- h. Familiar with the statistical and data analysis techniques, such as Matlab, ImageJ, Origin, Excel, and etc.
- i. Had some experience with other *in vivo* biomedical imaging techniques, such as the optical imaging and micro-CT imaging.

#### 4. POSTGRADUATE TRAINING

**Nov. 8<sup>th</sup>, 1997 – Mar. 31<sup>st</sup>, 2003** Postdoc Associate of Dr. A. Dean Sherry group, Department of Chemistry University of Texas at Dallas  
“Lanthanide complexes for biomedical usages”

#### 5. STAFF AND FACULTY POSITIONS

**Jun. 9<sup>th</sup>, 2019 – current** Senior Manager, Pre-Clinical Imaging Core, Center for Biometric Analysis, The Jackson Laboratory, Bar Harbor, ME, USA

**Dec. 9<sup>th</sup>, 2019 – Jun. 8<sup>th</sup>, 2021** Pre-Clinical Imaging Scientist, Center for Biometric Analysis, The Jackson Laboratory, Bar Harbor, ME, USA

**Jan. 1<sup>st</sup>, 2013 – Dec. 4<sup>th</sup>, 2019** Manager of UT-SAIR (Small Animal Imaging Resource), Sr. Research Engineer, Mouse MRI Core, Advanced Imaging Research Center, University of Texas Southwestern Medical Center, Dallas, TX, USA

**Jan. 3<sup>rd</sup>, 2011 – Dec. 30<sup>th</sup>, 2012** Research Scientist, Simons Comprehensive Cancer Center, University of Texas Southwestern Medical Center, Dallas, TX  
“New Responsive Contrast agents for MRI Cancer imaging”

**Sep. 1<sup>st</sup>, 2004 – Dec. 15<sup>th</sup>, 2010** Research Assistant Professor, Department of Radiology University of Washington, Seattle, WA, USA  
“Novel contrast agents for multi-modality applications”

**Apr. 1<sup>st</sup>, 2003 – Aug 27<sup>th</sup>, 2004** Research Scientist, Department of Radiology University of Texas Southwestern Medical Center, Dallas, TX  
“Novel contrast agents for MRI molecular imaging”

**Oct. 16<sup>th</sup>, 1995 – Oct. 30<sup>th</sup>, 1997** Assistant Researcher, Changchun Institute of Applied Chemistry, Chinese Academy of Sciences, P. R. China  
“NMR studies of lanthanide complexes”

#### 6. HONORS

**May, 2002** E. K. Zavoisky Stipend, International Society of Magnetic Resonance in Medicine

**Oct, 1995** Excellent Doctoral Fellow, Chinese Academy of Sciences

## 7. PROFESSIONAL ORGANIZATIONS

- Jan. 1<sup>st</sup>, 2021 – current** Member of International Society of Magnetic Resonance in Medicine (ISMRM)
- Jan. 1<sup>st</sup>, 2016 – Dec. 30<sup>th</sup>, 2016** Member of International Society of Magnetic Resonance in Medicine (ISMRM)
- Jan. 1<sup>st</sup>, 2004 – Dec. 30<sup>th</sup>, 2010** Member of International Society of Magnetic Resonance in Medicine (ISMRM)

## 8. RESEARCH FUNDING

### Accomplished:

- 02/01/2007 – 03/31/2011** PI, Department of Defense, US Army, “Dual-modality contrast agents for MRI and micro-CT imaging of prostate cancer”
- 01/01/2009 – 06/30/2010** PI, Provost Bridge Fund, UW, “Smart MRI Contrast Agent”
- 05/01/2006 – 04/30/2008** PI, Royalty Research Funds, UW, “Novel MRI Thermometry Based upon PARACEST Agents”

## 9. BIBLIOGRAPHY (*Selected among 83+ peer-reviewed articles*)

- (1) Selvaraj, UM, et al. *Brain, Behavior, and Immunity*, **2021**, 95: 502-13. “Delayed diapedesis of CD8 T cells contributes to long-term pathology after ischemic stroke in male mice” (I am the 6<sup>th</sup> coauthor).
- (2) Tiwari, V., et al. *Magn. Reson. Med.*, **2020**, 84(3): 1152-60. “In vivo MRS measurement of 2-hydroxyglutarate in patient-derived IDH-mutant xenograft mouse models versus glioma patients” (I am the 8<sup>th</sup> coauthor).
- (3) Ratnaker, S.J., et al. *Angew Chem. Int. Ed.*, **2020**, 59(48): 21671-76. “A Frequency-Selective pH-Responsive paraCEST Agent”. (I am the 4<sup>th</sup> coauthor)
- (4) Li, H-D., et al. *JCI Insight*, **2020**, 5(14): e138829. “A PoleP286R mouse model of endometrial cancer recapitulates high mutational burden and immunotherapy response”. (I am the 10<sup>th</sup> coauthor)
- (5) Gong, K. et al. *Nat. Cancer*, **2020**, 1: 394-409. “EGFR inhibition triggers and adaptive response by co-opting antiviral signaling pathways in lung cancer”. (I am the 14<sup>th</sup> coauthor).
- (6) Carter, A.M., et al. *PNAS*, **2020**, 117(31): 18401-11. “Phosphoprotein-based Biomarkers as Predictors for Cancer Therapy” (I am the 9<sup>th</sup> coauthor)
- (7) Kikuchi, K.; Ishimatsu, K.; Zhang, S.; Dimitrov, I.E.; Honda, H.; Sherry, A.D.; Takahashi, M. *Contrast Media & Molecular Imaging*, **2018**, 3141789. “Presaturation Power Adjusted Pulsed CEST: A Method to Increase Independence of Target CEST Signal”
- (8) Martins, A.F.; Jordan, V.C.; Bochner, F.; Chirayil, S.; Parawithana, N.; Zhang, S.; Lo, S-T.; Wen, X.; Zhao, P.; Neeman, M.; Sherry, A.D. *JACS*, **2018**, 140(50): 17456 - 64. “Imaging Insulin Secretion from Mouse Pancreas by MRI Is Improved by Use of a Zinc-Responseive MRI Sensor with Lower Affinity for Zn<sup>2+</sup> Ions”
- (9) Gong, K., et al. *J. Clin. Invest.*, **2018**, 128(6): 2500 - 18. “TNF-driven adaptive response mediates resistance to EGFR inhibition in lung cancer” (I am the 14<sup>th</sup> coauthor)

- (10) Selvaraj, U. M.; Zuurbier, K. R.; Whoolery, C. W.; Plautz, E. J.; Chambliss, K. L.; Kong, X.; Zhang, S.; Kim, S. H.; Katzenellenbogen, B. S.; Katzenellenbogen, J. A.; Mineo, C.; Shaul, P. W.; Stowe, A. M. *Endocrinology*, **2018**, 159(11): 3848 - 59. "Selective Nonnuclear Estrogen Receptor Activation Decreases Stroke Severity and Promotes Functional Recovery in Female Mice"
- (11) Park, S.K.; Rosenthal, T.R.; Williams, J.S.; Shelton, J.M.; Takahashi, M.; Zhang, S.; Bobulescu, I.A. *J. Invest. Med.*, **2018**, 66: 1037 - 44. "Metabolic and Cardiovascular Effects of Chronic Mild Hyperuricemia in Rodents"
- (12) Hajarnis, S.; Lakhia, R.; Yheskel, M.; Williams, D.; Sorourian, M.; Liu, X.; Aboudehen, K.; Zhang, S.; Kersjes, K.; Galasso, R.; Li, J.; Kaimal, V.; Lockton, S.; Davis, S.; Flaten, A.; Johnson, J. A.; Holland, W. L.; Kusminski, C. M.; Scherer, P. E.; Harris, P. C.; Trudel, M.; Wallace, D. P.; Igarashi, P.; Lee, E. C.; Androsavich, J. R.; Patel, V. *Nat. Commun.*, **2017**, 8: 14395. "MicroRNA-17 Family Promotes Polycystic Kidney Disease Progression through Modulation of Mitochondrial Metabolism"
- (13) Yu, M.; Zhou C.; Liu, L.; Zhang, S.; Sun, S.; Hankins, J.; Sun, X.; Zheng, J. *Angew. Chem. Int. Ed.*, **2017**, 56(15), 4314 - 9. "Interactions of Renal-Clearable Gold Nanoparticles with Tumor Microenvironments: Vasculature and Acidity Effects"
- (14) Moon, J.; Zhou, H.; Zhang, L.; Tan, W.; Liu, Y.; Zhang, S.; Morlock, L. K.; Bao, X.; Palecek, S. P.; Feng, J.; Williams, N. S.; Amatruda, J. F.; Olson, E. N.; Bassel-Duby, R.; Lum, L. *PNAS*, **2017**, 114(7): 1649 - 54. "Blockade to Pathological Remodeling of Infarcted Heart Tissue Using a Porcupine Antagonist"
- (15) Hussain R.Z.; Miller-Little, W.A.; Lambrach-Washington D.; Jaramillo, T. C.; Takahashi, M.; Zhang, S.; Fu, M.; Cutter, G. R.; Hayardeny, L.; Powell, C. M.; Rosenberg, R. N.; Stüve, O. *J. Neuroimmunol.* **2017**, 300: 100 - 9. "Laquinimod has no effects on brain volume or cellular CNS composition in the F1 3xTg-AD/C3H mouse model of Alzheimer's disease"
- (16) Zhang, L.; Martins, A. F.; Mai, Y.; Zhao, P.; Funk, A. M.; Clavijo Jordan, M. V.; Zhang, S.; Chen, W.; Wu, Y.; Sherry, A. D. *Chem Eur. J.*, **2017**, 23(8): 1752 - 8. "Imaging Extracellular Lactate In Vitro and In Vivo Using CEST MRI and a Paramagnetic Shift Reagent"
- (17) Singh, J.; Rustagi, V.; Zhang, S.; Sherry, A. D.; Udugamasooriya, D. G. *Magn. Res. Chem.* **2017**, 55(8): 747 - 53. "On-bead combinatorial synthesis and imaging of europium(III)-based paraCEST agents aids in identification of chemical features that enhance CEST sensitivity"
- (18) Clavijo Jordan, M. V.; Lo, S-T.; Chen, S.; Preihs, C.; Chirayil, S.; Zhang, S.; Kapur, P.; Li, W-H.; De Leon-Rodriguez, L. M.; Lubag, A. J. M.; Rofsky, N. M.; Sherry, A. D. *PNAS*, **2016**, 113(37): E5464-E5471. "Zinc-Sensitive MRI Contrast Agent Detects Differential Release of Zn(II) Ions from the Healthy vs. Malignant Mouse Prostate"
- (19) Preslar, A. T.; Tantakitti, F.; Park, K.; Zhang, S.; Stupp, S. I.; Meade, T. J. *ACS Nano*, **2016**, 10(8): 7376 - 84. "<sup>19</sup>F Magnetic Resonance Imaging Signals from Peptide Amphiphile Nanostructures Are Strongly Affected by Their Shape"
- (20) Wu, Y.; Zhang, S.; Soesbe, T. C.; Yu, J.; Vinogradov, E.; Lenkinski, R. E.; Sherry, A. D. *Magn. Reson. Med.*, **2016**, 75(6): 2432 - 41. "pH imaging of mouse kidneys in vivo using a frequency-dependent paraCEST agent"
- (21) Yousuf, M. A.; Tan, C.; Torres-Altora, M. I.; Lu, F-M.; Plautz, E.; Zhang, S.; Takahashi, M.; Hernandez, A.; Kernie, S. G.; Plattner, F.; Bibb, J. A. *J. Neurochem.*, **2016**, 138: 317 - 27. "Involvement of aberrant cyclin-dependent kinase 5/p25 activity in experimental traumatic brain injury"
- (22) Jiang, W.; Lumata, L.; Chen, W.; Zhang, S.; Kovacs, Z.; Sherry, A.D.; Khemtong C. *Sci. Reports*, **2015**, srep09104. "Hyperpolarized <sup>15</sup>N-pyridine Derivatives as pH-Sensitive MRI Agents"
- (23) Kumar, A.; Zhang, S.; Hao, G.; Hassan, G.; Ramezani, S.; Sagiya, K.; Lo, S-T.; Takahashi, M.; Sherry, A.D.; Oz, O.K.; Kovacs, Z.; Sun, X. *Bioconj. Chem.*, **2015**, 26(3): 549 - 58. "Molecular Platform for Design and Synthesis of Targeted Dual-Modality Imaging Probes"

- (24) Soesbe, T.C.; Ratnaker S.J.; Mile, M.; Zhang, S.; Do, Q.N.; Kovacs, Z.; Sherry A.D. *Magn. Res. Med.*, **2014**, 71(3): 1179 - 85. "Maximizing T<sub>2</sub>-exchange in Dy<sup>3+</sup>DOTA-(amide)<sub>x</sub> chelates: Fine-tuning the water molecule exchange rate for enhanced T<sub>2</sub> contrast in MRI"
- (25) Zhang, S.; Zhou, K.; Huang, G.; Takahashi, M.; Sherry, A. D.; Gao, J. *Chem. Commun.*, 2013, 49(57): 6418-6420. "A novel class of polymeric pH-responsive MRI CEST agents"
- (26) Huang, X.; Huang, G.; Zhang, S.; Sagiya, K.; Togao, O.; Ma, X.; Wang, Y.; Li, Y.; Soesbe, T.C.; Sumer, B.; Takahashi, M.; Sherry, A.D.; Gao, J. *Angew Chem. Int. Ed.*, **2013**, 52(31): 8074 - 78. "Multi-Chromatic pH-Activatable F-19-MRI Nanoprobes with Binary ON/OFF pH Transitions and Chemical-Shift Barcodes"
- (27) Zhou, K.; Liu, H.; Zhang, S.; Huang, X.; Wang, Y.; Huang, G.; Sumer, B.; Gao, J. *Journal of the American Chemical Society*, **2012**, 134(18): 7803 - 11. "Multicolored pH-Tunable and Activatable Fluorescence Nanoplatfrom Responsive to Physiologic pH Stimuli"
- (28) Bui, T.; Stevenson, J.; Hoekman, J.; Zhang, S.; Maravilla, K.; Ho, R.J.Y. *PLoS One*, **2010**, 5(9): e13082. "Novel Gd Nanoparticles Enhance Vascular Contrast for High-Resolution Magnetic Resonance Imaging"
- (29) Ren, J.; Trokowski, R.; Zhang, S.; Malloy, C. R.; Sherry, A. D. *Magnetic Resonance in Medicine* **2008**, 60, 1047 - 55. "Imaging the tissue distribution of glucose in livers using a PARACEST sensor"
- (30) Park, I. K., Ng, C. P., Wang, J., Chu, B., Yuan C., Zhang, S., Pun, S.H. *Biomaterials*, **2008**, 29: 724 - 32. "Determination of nanoparticle vehicle unpackaging by MR imaging of a T<sub>2</sub> magnetic relaxation switch"
- (31) Woods, M.; Caravan, P.; Gerald, C.F.G.C.; Greenfield, MT.; Kiefer, GE.; Lin, M.; McMillan, K.; Prata, MIM.; Santos, AC.; Sun, X.; Wang, J.; Zhang, S.; Zhao, P.; and Sherry, AD. *Investigative Radiology*, **2008**, 43(12): 861 - 70, "The Effect of the Amide Substituent on the Biodistribution and Tolerance of Lanthanide(III) DOTA-Tetraamide Derivatives"
- (32) Vasalatiy, O.; Zhao, P.; Zhang, S.; Aime, S.; Sherry, AD.; *Contrast Media & Molecular Imaging*, **2006**, 1: 10-14. "Catalytic Effects of Apoferritin Interior Surface Residues on Water Proton Exchange in Lanthanide Complexes"
- (33) Garcia-Martin, ML.; Martinez, GV.; Raghunand, N.; Sherry AD.; Zhang, S.; Gillies, RJ.; *Magnetic Resonance in Medicine*, **2006**, 55(2): 309 - 15. "High Resolution pH<sub>e</sub> Imaging of Rat Glioma Using pH-Dependent Relaxivity"
- (34) Zhang, S., Malloy, CR., Sherry, AD. *Journal of the American Chemistry Society*, **2005**, 127: 17572 - 3. "MRI Thermometry Based on PARACEST Agents"
- (35) Zhang, S., Jiang X, Sherry, AD. *Helvetica Chimica Acta*, **2005**, 88: 923 - 35. "Modulation of the lifetime of water bound to lanthanide metal ions in complexes with ligands derived from 1,4,7,10-tetraazacyclododecane tetraacetate (DOTA)"
- (36) Sherry, AD, Zhang, S., Woods, M. *ACS Symposium Series – Medical Inorganic Chemistry*, **2005**, 903: 151 - 65. "Water Exchange Is the Key Parameter in the Design of Next-generation MRI Agents"
- (37) Vinogradov, E., Zhang, S., Lubag, A., Balschi, J.A., Sherry, A.D., Lenkinski, R.E. *Journal of Magnetic Resonance*, **2005**, 176: 54-63. "On-resonance low B<sub>1</sub> pulses for imaging of the effects of PARACEST agents"
- (38) Woessner DE, Zhang S., Merritt ME, Sherry, AD. *Magnetic Resonance in Medicine* **2005**, 53(4): 790 - 9. "Numerical solution of the Bloch equations provides insights into the optimum design of PARACEST agents for MRI"
- (39) Woods, M.; Zhang, S.; Sherry, A.D. *Current Medicinal Chemistry: Immunology, Endocrine & Metabolic Agents*, **2004**, 4(4), 349 - 69. "Toward the design of MR agents for imaging β-cell function"

- (40) Trokowski, R.; Zhang, S.; Sherry, A. D. *Bioconjugate Chemistry*, **2004**, 15(6), 1431- 40. "Cyclen-Based Phenylboronate Ligands and Their Eu<sup>3+</sup> Complexes for Sensing Glucose by MRI"
- (41) Zhang, S.; Trowkoski, R.; Sherry, A.D. *Journal of the American Chemistry Society*, **2003**, 125, 15288 - 9. "A paramagnetic CEST agent for imaging glucose by MRI"
- (42) Zhang, S.; Merritt, M.; Weossner, D.E.; Lenkinski, R.E.; Sherry, A.D. *Accounts of Chemical Research*, **2003**, 36, 783-790. "PARACEST Agents: Modulating MRI Contrast via Water Proton Exchange"
- (43) Zhang, S.; Sherry, A.D. *J. Solid State Chem.*, **2003**, 171(1-2): 38-43. "Physical Characteristics of Lanthanide Complexes that act as Magnetization Transfer (MT) Contrast Agents"
- (44) Woods, M.; Zhang, S.; Ebron, V.H.; Sherry, A.D. *Chem. Eur. J.*, **2003**, 9(19), 4634-40. "pH-sensitive modulation of the second hydration sphere in lanthanide (III) tetraamide DOTA complexes: a novel approach to smart MR contrast media"
- (45) Woods, M.; Kovacs, Z.; Zhang, S.; Sherry, A.D. *Angewandte Chemistry International Edition*, **2003**, 42(47), 5889-5892. "Towards the Rational Design of MRI Contrast Agents; Isolation of the Two Coordination Isomers of Lanthanide DOTA-type Complexes"
- (46) Raghunand, H.; Howison, C.; Sherry, A.D.; Zhang, S.; Gillies, R.J. *Magnetic Resonance in Medicine*, **2003**, 49(2): 249-257. "Renal and Systemic pH Imaging by Contrast-enhanced MRI"
- (47) Geraldes, C.G.F.C.; Zhang, S.; Sherry, A.D. *Bioinorganic Chemistry and Applications*, **2003**, 1(1): 1-23. "Comparison of Crystal Field Dependent and Independent Methods to Analysis Lanthanide Induced Shifts in Axially Symmetric Complexes: Part II. Systems with a C<sub>4</sub> symmetry"
- (48) Alves, F.C.; Donato, P.; Sherry, A.D.; Zaheer, A.; Zhang, S.; Lubag, A.J.M.; Merritt, M.E.; Lenkinski, R.E.; Frangioni, J.V.; Neves, M.; Prata, M.I.M.; Santos, A.C.; De Lima, J.J.P.; Geraldes, C.G.F.C. *Investigative Radiology*, **2003**, 38(12): 750-60. "Silencing of phosphonate-gadolinium magnetic resonance imaging contrast by hydroxyapatite binding"
- (49) Zhang, S.; Michaudet, L.; Burgess, S.; Sherry, A.D. *Angewandte Chemistry International Edition*, **2002**, 41(11): 1919-1921. "The amide protons of a ytterbium(III)-DOTA-tetraamide complex act as efficient antennae for transfer of magnetization to bulk water"
- (50) Zhang, S.; Wu, K.; Sherry, A.D. *Journal of the American Chemistry Society*, **2002**, 124(16): 4226-4227. "Unusually Sharp Dependence of Water Exchange Rate versus Lanthanide Ionic Radii for a Series of Tetra-amide Complexes"
- (51) De Leon, L.M.; Ortiz, A.; Weiner, A.L.; Zhang, S.; Kovacs, Z.; Kodadek, T.; Sherry, A.D. *Journal of the American Chemistry Society*, **2002**, 124(14): 3514-3515. "Magnetic Resonance Imaging Detects a Specific Peptide-protein Binding Event"
- (52) Vander Elst, L.; Zhang, S.; Sherry, A.D.; Laurent, S.; Botterman, F.; Muller, R.N. *Academica Radiology*, **2002**, 9: s297-s299. "Dy-complexes as High Field T<sub>2</sub> Contrast Agents: Influence of Water Exchange Rates"
- (53) Raghunand, N.; Zhang, S.; Sherry, A.D.; Gilles, R.J. *Academica Radiology*, **2002**, 9: s481-s483. "In Vivo Magnetic Resonance Imaging of Tissue pH by Using a Novel pH-Sensitive Contrast Agent, GdDOTA-4AmP"
- (54) Ren, J.; Zhang, S.; Sherry, A.D.; Geraldes, C.F.G.C. *Inorganic Chimica Acta*, **2002**, 339: 273-282. "Analysis of lanthanide induced NMR shifts using a crystal field independent method: application to complexes of tetraazamacrocyclic ligands"

- (55) Zhang, S.; Winter, P.; Wu, K.; Sherry, A.D. *Journal of the American Chemistry Society*, **2001**, 123(7): 1517-1518. "A Novel Europium-based MRI Contrast Agent"
- (56) Zhang, S.; Wu, K.; Sherry, A.D. *Investigative Radiology*, **2001**, 36: 82-86. "Gd<sup>3+</sup> Complexes with Slowly Exchanging Bound Water Molecules May Offer Advantages in the Design of Functional MRI Agents"
- (57) Zhang, S.; Wu, K.; Beaver, M.; Sherry, A.D. *Inorganic Chemistry*, **2001**, 40: 4284-4291. "<sup>1</sup>H and <sup>17</sup>O NMR Detection of A Lanthanide-bound Water Molecule at Ambient Temperatures in Pure Water as Solvent"
- (58) Zhang, S.; Kovacs, Z.; Burgerss, S.; et al. *Chemistry -- A European Journal*, **2001**, 7: 288-296. "Lanthanide-DOTA-Bis(amide) Complexes: Direct NMR Evidence for Differences in Water Molecule Exchange Rates for Coordination Isomers"
- (59) Li, X.; Zhang, S.; Zhao, P.; Kovacs, Z.; Sherry, A.D. *Inorganic Chemistry*, **2001**, 40: 6572-6579. "Synthesis and NMR Studies of New DOTP-like Lanthanide (III) Complexes Containing a Hydrophobic Substituent on One Phosphonate Side Arm"
- (60) Zhang, S.; Wu, k.; Sherry, A.D. *Angewandte Chemistry International Edition*, **1999**, 38(21): 3192-3194. "A novel pH-sensitive MRI contrast agent"
- (61) Zhang, S.; Sun, D.; Li, X.; Pei, F.; Liu, S. *Fullerene Science and Technology*, **1997**, 5(7): 1635-1643. "Synthesis and solvent enhanced relaxation property of water-soluble endohedral metallofullerenols"

#### 10. Others (Patents):

- (1) Sherry, A.D; Zhang, S.; Wu, K. *US Patent*, US6,746,662 (Date of Patent: Jun. 8, **2004** "pH-sensitive MRI contrast agents"
- (2) Sherry, A.D; Zhang, S.; Wu, K. *US Patent*, US6,875,419 (Date of Patent: Apr. 5, **2005** "Paramagnetic metal ion-based macrocyclic magnetization transfer contrast agents and method of use"